

# Acuson Maple



- Console System
- Color, Power, PW, Triplex
  - Calculation Packages
- Speckle Reduction, Harmonics
  - Adjustable Height
  - Touch Screen
- HDMI, S video, Networkable
  - Backlit Keyboard
- 21.5 in LED Monitor

## Specifications

Scanning Method	Convex, Linear, Endocavity, Phased,
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power, 2D/CFM
Focal Zones	up to 8 Adjustable
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Speckle Reduction TCE, Edge Enhancement, Smoothing, Zoom, Compounding
Image Filing	Hard Drive 500GB, 6 USB, CD/DVD, Dicom Option, Networking
Measurements	Distance, Circ., Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, Card, OB/GYN, Uro, Thyroid, Venous, TCD, Emergency
Monitor	21.5 inch 1080 Flat screen LED
Probe Connections	3, 4 optional
Dimensions	H 55.2-61.6 x W 23.1 x D 31.2"
Weight	Approx. 144 Lbs.

## Probes and Applications

Convex 5C1a 5C-2v 7C2 11M3	1.4-5.0 MHz. 1.7-4.9 MHz 2.1-7.0 MHz 3.5-11 MHz	OBGYN, Abd, Card, Abd, MSK Small Parts
Endocavity 9VEV 10MC3	3.1-8.7 MHz 3.5-10.2MHz	Abd, OB, Card
4D 9VC2	1.8-8.5 MHz	OB 4D
Linear 14L4a 16L4 L10-5v	4-12.7 MHz 4.3-15.4MHz 4.2 -12 MHz	Breast, Thy, SP MSK, Vasc, Uro
Phased CW2 CW5	2 MHz 5 mhZ	Card Vasc

### Options:

Gel Warmer, Pullout Keyboard, Battery, 3D/4D, Auto OB, Auto Fol, Vascular, Auto OB, CIMT, Stress, Biopsy Guides, Elast, Siescan, Printers, DICOM, Needle, Modality Worklist, Structured Reporting



## DIAGNOSTIC HEALTH GROUP

222 Institute Street • Smithfield, VA 23430 •

E-Mail: [waxmed@verizon.net](mailto:waxmed@verizon.net)

800 669-3442 • Fax 757 357-585 • [www.dhgsommed.com](http://www.dhgsommed.com)