

Acuson NX3



- 10.4 in Touch Screen
- Color, Power, PW Triplex
 - Calculation Packages
- Speckle Reduction, Harmonics
 - Auto Trace & Auto Calc
- Can use Legacy Transducers
 - VGA, DVI
 - Backlit Keyboard
 - 21.5 in Monitor

Specifications

Scanning Method	Convex, Linear, Endocavity, Phased,
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power, 2D/CFM
Focal Zones	up to 4 Adjustable span & depth
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Speckle Reduction TCE, Edge Enhancement, Smoothing, Zoom, Compounding
Image Filing	Hard Drive, 6 USB, CD/DVD, Dicom Option, Networking
Measurements	Distance, Circ., Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, Card, OB/GYN, Uro, Thyroid, Venous, TCD, Emergency
Monitor	21.5 inch 1080 Flat screen
Probe Connections	3, 4 optional
Dimensions	H 58.9 x W 21.3 x D 36.3 inches
Weight	Approx. 165 Lbs.

Probes and Applications

Linear VF10-5 VF12-4 VF13-5 VF16-5	5-10 MHz. 4-12 MHz 5-13 MHz 5-16 MHz	Breast, Thyroid, MSK, Vasc, Small Parts
Convex CH5-2 C8-5	2.5,3.6,5MHz 3.1,3.6,4.4,5 3.1-8.8 MHz	Abd, OB, Card
Endocavity EC9-4 EC10-5	5,6.7,7.3MHz 5.7,6.6,9.4 5.7-10 MHz	GYN, Urology
Phased P4-2 P8-4	2,-3.6 MHz 4-8 MHz	Cardiac, TCD Ped Card
4D C8F8	3.3,4.2,5MHz 4.4,5,5.7MHz	OB

Options:

HD Zoom, Quick Start, 3D/4D, Vascular Enhancement, Auto OB, CIMT, Stress, Biopsy Guides, Printers, DICOM, Modality Worklist, Structured Reporting, Gel Warmer, Adjustable Keyboard



DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: DHGsonomed@aol.com
800 669-3442 • Fax 301 695-7940 • www.dhgsonomed.com