

Apogee 1000



- Laptop style
- Color, Power, PW Triplex
 - Calculation Packages
- Speckle Reduction, Harmonics
 - Auto Trace & Auto Calc
- 5 Frequency Transducers
 - Backlit Keyboard
- 15" Adjustable Display

Specifications

Scanning Method	Convex, Linear, Endocavity, Phased
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power
Focal Zones	up to 8 Adjustable span & depth
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Speckle Reduction, Edge Enhancement, Smoothing, Zoom
Image Filing	500 GB Hard Drive, BMP, JPG USB
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid
Monitor	15 inch LCD Adjustable
Probe Connections	1
Dimensions	H 13.94 x W 19.37 x D 8.3 inches
Weight	Approx. 13.1 Lbs.

Probes and Applications

Linear L8LC L10LC L5LC	MHz 5,6,6,7.5,10,12 7,9,10,12,14 4,4,7,5.5,6,7,2	25,38&50 mm Breast, Thyroid, MSK, Vasc, Small Parts
Convex C3LC R60 C6LC 17.5 C3IC R20	MHz3.3,4.2,5 4.5,5.2,6,6.7,7.5 2,2.5,3.3,4.2,5	Abd, OB, Card Ped Abd, Ped Car Card
Endocav V6LC U5LC ECBP	MHz 4,5,6,7,9 MHz 4,4,7,5.5,6.2,7 5,6,6,7.5,10,12	GYN, Urology Uro Linear Uro BiPlane
4D 4DL40J V3120J	MHz 4,4,7,5.5,6.2,7 4,5,6,7,9	OB Endovag
Phased P3FC P5FC	MHz 1.7,2,2.5,3.3,4 4,4,7,5.5,6.2,7	Card,Abd,TCD Ped Card,Abd

Options:

Printer, Footswitch, Dicom,
Biopsy Guides, Elastography, WIFI,
Battery, Cart, MultiProbe 2 or 4



DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net
800 669-3442 • Fax 757 357-5985 • www.dhgsommed.com