

# Apogee 1200 T



- Portable Lunchbox style
- Color, Power, PW Triplex
  - Calculation Packages
- Speckle Reduction, Harmonics
  - Auto Trace & Auto Calc
  - 5 Frequency Transducers
- VGA, S video, HDMI outputs
  - Backlit Keyboard
  - Upgradeable to 4D

## Specifications

Scanning Method	Convex, Linear, Endocavity, BiPlane,
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power
Focal Zones	up to 8 Adjustable span & depth
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Speckle Reduction, Edge Enhancement, Smoothing, Zoom
Image Filing	Hard Drive, USB, Dicom Option
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid
Monitor	15 inch LCD
Probe Connections	2
Dimensions	H 14.5 x W 15.3 x D 8.2 inches
Weight	Approx. 24 Lbs.

## Probes and Applications

Linear L8L38J L8L50J L10L38J	5,6,6,7.5,10, 12 MHz. 8,9,10,11,12	Breast, Thyroid, MSK, Vasc, Small Parts
Convex C3L60J C3L20J C5L20J	2,2.5,3.3,4.2 ,5 MHz 4,4.7,5.5,6.2 ,7 MHz.	Abd, OB, Card  Ped Abd, Ped Card
Endocavity V6L11J U5L50J	4,5,6,7,9 MHz. 5,6,6,7.5,10 12 MHz.	GYN, Urology
4D 4DL40J V3I20J	4,4.7,5.5,6.2- ,7 MHz 4,5,6,7,9MHz	OB Endovag

### Options:

Printer, Footswitch, Dicom,  
Biopsy Guides, Elastography, WIFI,  
3D/4D



## DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: [DHGsonomed@aol.com](mailto:DHGsonomed@aol.com)  
800 669-3442 • Fax 301 695-7940 • [www.dhgsonomed.com](http://www.dhgsonomed.com)