

SIUI Apogee 2100



- Portable Lunchbox style
- Angle Adjustable Display
- Color, Power, PW Triplex
 - Calculation Packages
- Speckle Reduction, Harmonics
 - 2 Active Probe Ports
 - 5 Frequency Transducers
- VGA, S video, HDMI outputs
 - Upgradeable to 4D

Specifications

Scanning Method	Convex, Linear, Endocavity, BiPlane,
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power
Focal Zones	up to 8 Adjustable span & depth
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Speckle Reduction, Edge Enhancement, Smoothing, Zoom
Image Filing	Hard Drive, USB, Dicom Option
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid
Monitor	15 inch LCD Angle Adjustable
Probe Connections	2
Dimensions	H 24.8 x W 14.17 x D 5.7 inches
Weight	Approx. 32 Lbs.

Probes and Applications

Linear L8LN 38or50mm L10LN 38	5,6,6,7.5,10, 12 MHz. 8,9,10,12,14	Breast, Thyroid, MSK, Vasc, Small Parts
Convex C3LN C6LN R17	2,2.5,3.3,4.2 ,5 MHz 4.5,5.2,6.0, 6.7,7.5MHz.	Abd, OB, Card Ped Abd, Ped Card
Endocavity V6LN U5LN	4,5,6,7,9 MHz. 4,4.7.5.5,6.2 7 MHz.	GYN, Urology Urology
P3FN	1.7,2,2.5,3.3, 5 MHz	Cardiac Abdominal TCD
4D C5LN C5LN Low	4,4.7,5.5,6.2 7 MHz 2,2.5,3.3,4.2 5MHz	OB

Options:

Printer, Footswitch, Dicom,
Biopsy Guides, Elastography, WIFI,
Cart, ECG, Case, 3D/4D, 4D Pro



DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net
800 669-3442 • Fax 757 357-5985 • www.dhgsomed.com