

# CTS 8800 Plus



- Portable Lunchbox style
- B&W, M, PW Doppler
- Calculation Packages
- Speckle Reduction, Trapezoidal
- Harmonics, Auto Adjust
- 5 Frequency Transducers
- VGA, S video & HDMI outputs
- Backlit Keyboard
- 250 GB HD, USB

## Specifications

Scanning Method	Convex, Linear, Endocavity, BiPlane, 4D
Display Modes	B, 2B, 4B, B/M, M, PW
Focal Zones	Adjustable
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Edge Enhancement, Smoothing, Zoom
Image Filing	Hard Drive, USB, Dicom Option
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, SD, OB/GYN, Uro
Monitor	15 inch LCD
Probe Connections	2
Dimensions	H 15.5 x W 14.6 x D 8.2 inches
Weight	Approx. 24 Lbs.

## Probes and Applications

Linear L7L38K-G L7150K-G L10L38K-G L10L25K-G	5-12 MHz. 5-12 MHz. 7-14 MHz. 7-14 MHz.	Breast, Thyroid, MSK, Vasc, Small Parts
Convex C3L60K-G C3L40K-G	2-5 MHz. 2-5 MHz.	Abd, OB, Card Ped Abd, Ped Card
Micro Conv C3I20K-G C6I15K-G	2-5 MHz. 4.5-7.5 MHz.	Card, Abd, Ped Card, Ped ABd
Endocavity V6L11K-G L5L50K-G	4,5,6,7,9 MHz. 4,4.7,5.5,6.2,7	GYN, Urology
4D 4DL40BK-G	2-5 MHz.	OB

Options:  
Printer, Footswitch, Dicom,  
Biopsy Guides, Compound Imaging,  
Panoramic, CFM, 4D, WIFI



## DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: [DHGsonomed@aol.com](mailto:DHGsonomed@aol.com)  
800 669-3442 • Fax 301 695-7940 • [www.dhgonomed.com](http://www.dhgonomed.com)