

# Mindray DP50 Expert



- Portable and Compact
- 15 inch LCD Monitor
- Programable Application Settings
- Multi-Frequency Transducers
  - Calculation Packages
  - Hard Drive and USB
- Speckle Reduction, Compounding
  - Harmonic, Auto Adjust
  - Battery
- Two Probe Connections

## Specifications

Scanning Method	Linear, Convex, Endocavity
Display Modes	B, Dual B, Quad B, B/M, M
Focal Zones	Adjustable
Image Adjustments	Gain, TGC 8 slides, Depth, Dynamic Range, Zoom, Edge Enhancement, Post Processing , Persistence
Image Filing	Hard Drive, 2 USB,
Measurements	Distance, Area, Circumference, Volume, OB/GYN, Urology, Abd, Small Parts
Monitor	15 inch with tilt
Probe Connections	2

## Probes and Applications

Convex 35C50EA	2-5.2 MHz.	Abd, OB/GYN
MicroConvex 65C15EA	4.2-9.5MHz.	Ped Abd
Linear 75L38EA 75L53EA 10L24EA	3.3-13 MHz. 4.5-10 MHz 5-10 MHz	Thyroid, MSK, Small Parts, Breast
Endocavity 65EC10EA 65EB10EA	4.3-10 MHz. 4.3-10 MHz.	Transvaginal Transrectal BP

### Options:

Video Printer, Footswitch, Cart,  
Biopsy Guides, Dicom, CFM, PW



## DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: [waxmed@verizon.net](mailto:waxmed@verizon.net)  
800 669-3442 • Fax 757 357-5985 • [www.dhgsomed.com](http://www.dhgsomed.com)