

Mindray DC40



- Adjustable Keyboard
- 13.3 in Touchscreen
- Color, Power, PW Triplex
 - Battery
- Calculation Packages
- Speckle Reduction
 - Compounding
 - Harmonics

Specifications

| | |
|-------------------|--|
| Scanning Method | Convex, Linear, Phased, Endocavity, 4D |
| Display Modes | B, 2B, 4B, B/M, M, PW, CFM, Power, Duplex, Triplex |
| Focal Zones | up to 4, 16 Adjustable |
| Image Adjustments | Gain, Depth, Dynamic Range, frequency, Focus, Persistence, GM, Tint, Speckle, Edge, Smooth, Sweep, Zoom |
| Image Filing | 1 TB Hard Drive, 7 USB, |
| Measurements | Distance, Circumference, Area, Ellipse, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid, Vasc |
| Monitor | 21.5 LED, 13.3 Touch Screen |
| Probe Connections | 4 |
| Dimensions | D 34.5 x W21.4 x H 44.6-59 inches |
| Weight | Approx. 75 Kgs. 166 lbs Main Unit |

Probes and Applications

| | | |
|---|--|---|
| Linear 7L4B 38 L14-6NE 38 L12-3E 38 L7-3 38 14-6 26 | MHz 3-14 4gen2har 3.3-16 4-2 3-13.5 4-2 2.7-8 4-2 3.5-16 4-2 | Pediatric, Small Parts, MSK, Vascular, Ortho, Nerve |
| Convex 3C5A 6C2 C6-2 D7-2E 4D D7-2 4D | 1.3-6 4-2 2.6-13.2 4-2 1.3-5.7 4-2 2-8.4 4-2 2-8.4 4-2 | OB/GYN, Ped, Abdominal, Vasc, Uro, Card. 4D |
| Endocavity V11-3 V10-4B V10-4 CB10-4E | 42.4-12.8 4-2 2.6-13.2 4-2 2.6-13.2 4-2 2.6-11 4-2 | OB/GYN, Uro BiPlane |
| Phased P4-2 P7-3 P10-4E | 1-5.2 4-2 2-8.2 4-2 3-11.8 4-2 | Card, Abd TCD, Ped |

Notes: Svideo, BNC, VGA, DVI outputs
Networkable, 8 TGC
Options: Printer, Footswitch, Physio,
Dicom, Biopsy Guides, IMT, Elast



DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net
800 669-3442 • Fax 757 357-5985 • www.dhgonomed.com