

# Mindray DC70



- Adjustable Keyboard
- 21.5 LED Monitor
- 13.3 in Touchscreen
- Color, Power, PW Triplex
  - Wireless
- Calculation Packages
- Speckle Reduction
  - Compounding
  - Harmonics

## Specifications

Scanning Method	Convex, Linear, Phased, Endocavity, 4D
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power, Duplex, Triplex
Focal Zones	up to 4, 16 Adjustable
Image Adjustments	Gain, Depth, Dynamic Range, Frequency, Focus, Persistence, GM, Tint, Speckle, Edge, Smooth, Sweep, Zoom
Image Filing	1 TB Hard Drive, 5 USB,
Measurements	Distance, Circumference, Area, Ellipse, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid, Vasc
Monitor	21.5 LED, 13.3 Touch Screen
Probe Connections	4
Dimensions	W21.7 x D 40.2 x H39.4 inches
Weight	Approx. 105 Kgs. 232 lbs Main unit

## Probes and Applications

Linear	MHz	
L14-6NE 38	5.4-16gen4har2	Pediatric, Small Parts, MSK, Vascular, Ortho, Nerve
L12-3E 38	3.3-16 4-2	
L9-3E 43	1.8-9.8 4-2	
L14-6WE51	2.7-8 4-2	
L14-3SWE54	3.5-16 4-2	
Convex		
SC6-1 60R	1.3-5.7 4-2	OB/GYN, Ped, Abdominal, Vasc, Uro, Card.
C5-1E1 60R	2.6-13.2 4-2	
C7-3 50R	2.6-7.2 4-2	
C11-3 15R	2.6-12.8 4-2	
SD8-1 45R	2.6-8.2 4-2	4D
D7-2E 40R	2.6-8.2 4-2	4D
Endocavity		
V11-3HE	2.6-12.8 4-2	OB/GYN, Uro
V11-4WE	2.6-12.8 4-2	
6LB7E	2.6-12.8 4-2	
CB10-4	2.6-12.8 4-2	
D11-3E	2.6-12.8 4-2	4D
Phased		
SP5-1E	1-5.2 3-3	Card, Abd, TCD
P4-2	1.3-4.6 3-3	
P7-3	2.3-7.2 3-3	Ped
P10-4E	3-11.4 3-3	Abd, Card
P7-3TE	2.3-7.2 3-3	Card, Ped

### Options:

Printer, Footswitch, Physio, Dicom, Biopsy Guides, IMT, Elast, Case, Smart Measures, Gel Warmer, Bar Codes, 4D



## DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net  
800 669-3442 • Fax 757 357-5985 • www.dhgsomed.com