

# Mindray M7



- Laptop style
- 15" LCD Display
- Color, Power, PW Triplex
  - Calculation Packages
  - iClear Speckle Reduction
  - iBeam compounding
  - iScape Panoramic
  - Harmonics
- Fast Power up, Simple Operation

## Specifications

Scanning Method	Convex, Linear, Endo, Phased, 4D
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power
Focal Zones	up to 16
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Edge Enhancement, Smoothing, Zoom
Image Filing	320GB Hard Drive, 2 USB, BMP, JPG, FRM, AVI, Dicom
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid
Monitor	15 inch LCD
Probe Connections	1 Multiprobe available on Cart
Dimensions	W 14.21 x D 14.06 x H 2.95 inches
Weight	Approx. 6.5 Kgs.

## Probes and Applications

Linear L14-6s L14-6N L12-4s L74s L7-3s 7L5s		Pediatric, Small Parts, MSK, Vascular, Ortho, Nerve
Convex C5-2s 6C2S 4CD4s	2-5 MHz. Microconvex 4D	OB/GYN, Ped, Abdominal, Vasc, Uro, Card.
Endocavity V10-4s V10-4Bs		OB/GYN, Uro
Phased P4-2s P7-3s P12-4		Card, Abd Ped TCD

Options:  
Printer, Footswitch, Cart, 3Port Ext,  
Biopsy Guides, IMT, WIFI,  
External Monitor,



## DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net  
800 669-3442 • Fax 757 357-5985 • www.dhgsomed.com