

# Mindray Z6



- Portable Lunchbox style
- 15" LCD Adjustable Display
- Color, Power, PW Triplex
  - Calculation Packages
  - 2 Active Probe Ports
  - iBeam compounding
    - iScape Panoramic
    - Harmonics
- VGA, S video, Compound Video

## Specifications

Scanning Method	Convex, Linear, Endocavity, Phased,
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power
Focal Zones	8
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Edge Enhancement, Smoothing, Zoom
Image Filing	Hard Drive, 4 USB, Dicom Option
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid
Monitor	15 inch LCD Angle Adjustable
Probe Connections	2
Dimensions	H 14.88 x W 16.34 x D 7.48 inches
Weight	Approx. 8.8 Kgs.

## Probes and Applications

Linear 7L4P 7L5P L14-6P 7LT4P	5-10 MHz 5-10 MHz 8-14 MHz 5-10 MHz	Pediatric, Small Parts, MSK, Vascular, Ortho, Nerve Interoperative
Convex CCSP 6C2P 6CV1P	2-5 MHz. 5-8.5 Micro 5-8.5 MHz	OB/GYN, Ped, Abdominal, Vasc, Uro, Card.
Endocavity V10-4BP CB10-4P 6LE7P	5-8.5 MHz 5-8.5 MHz 5-8.5 MHz	GYN Urology Linear Endo
Phased 2P2P	2-3.5 MHz	Card, Abd

Options:  
Printer, Footswitch, Dicom,  
Biopsy Guides, Elastography, WIFI,  
Cart, ECG, Case



## DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net  
800 669-3442 • Fax 757 357-5985 • www.dhgsomed.com