

Mindray Z60



- Portable Lunchbox style
- 15" LCD Adjustable Display
- Color, Power, PW Triplex
 - Calculation Packages
 - 2 Active Probe Ports
 - iBeam compounding
 - iScape Panoramic
 - Harmonics
- VGA, S video, Compound Video

Specifications

Scanning Method	Convex, Linear, Endocavity, Phased,
Display Modes	B, 2B, 4B, B/M, M, PW, CFM, Power
Focal Zones	8
Image Adjustments	Gain, Depth, Dynamic Range, Persistence, Line Density, Edge Enhancement, Smoothing, Zoom
Image Filing	50GB Hard Drive, 4 USB, Dicom Option
Measurements	Distance, Circumference, Area, Volume, Time, HR, Velocity Acc, RI, PI, SD, OB/GYN, Uro, Thyroid
Monitor	15 inch LCD Angle Adjustable
Probe Connections	2
Dimensions	H 14.1 x W 14.3 x D 5.8 inches
Weight	Approx. 8.2 Kgs.w/o Battery

Probes and Applications

Linear 7L4BP 7L5P L14-6P 7LT4P	5-10 MHz 5-10 MHz 8-14 MHz 5-10 MHz	Pediatric, Small Parts, MSK, Vascular, Ortho, Nerve Interoperative
Convex C6-2P 6C2P 6CV1P D6-2P	2-5 MHz. 5-8.5 Micro 5-8.5 MHz 2-5 MHz	OB/GYN, Ped, Abdominal, Vasc, Uro, Card 4D
Endocavity V10-4BP CB10-4P 6LE7P	5-8.5 MHz 5-8.5 MHz 5-8.5 MHz	GYN Urology Bi-Plane Linear Endo
Phased 2P2P P7-3P	2-3.5 MHz	Card, Abd Pediatric

Options:

Printer, Footswitch, Dicom, Biopsy Guides, Elastography, WIFI, Cart, ECG, Case, 4D



DIAGNOSTIC HEALTH GROUP

222 Institute St • Smithfield, VA 23430 • E-Mail: waxmed@verizon.net
800 669-3442 • Fax 757 357-5985 • www.dhgsomed.com