

Siemens Acuson P300



- Portable
- 15 inch XVGA Flat Screen
- Xview Speckle Reduction
- Hard Drive, CD, and USB Storage
- Networkable
- Console Performance
- Loaded with Features
- CFM, Power and PW Doppler

Specifications

Scanning Method	Linear, Convex, Phased Array Pencil Doppler
Display Modes	B, Dual B, B/M, M, PWD, Color, Power
Focal Zones	Up to 8
Image Adjustments	Image Gain, Doppler Gain, 7 TGC, Dynamic Range, Frequency, Gray Map, Colorize, Sharpness, Persistence, Harmonics
Image Filing	HD, CD, USB, Image, Clips (option)
Measurements	Distance, Ratio, Length, Ellipse, Area, Volume
Monitor	15 inch XVGA Flat Screen
Probe Connections	2 Imaging, 1 Doppler Probe
Dimensions	Closed 14 x 6.2 x 19.3 in. Open 14 x 16 x 19.3 in.
Weight	20 lbs. plus optional cart

Probes and Applications

Linear	10-18 MHz 6.6-13 MHz 3.3-7.5 MHz 5-12 MHz	MSK, Breast, Vascular, Small Parts
Convex	1.8 - 5 MHz 5 - 7.5 MHz	OB/GYN, Abd, Ped, Vasc.
Endocavity	5 - 7.5 MHz	Transvaginal, Transrectal
Pencil	5 MHz	Vascular

Options:

B&W Video Printer, Color Video Printer,
Inkjet Printer, VCR, Network, Dicom,
Carts, Biopsy Guides



DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: DHGsonomed@aol.com
800 669-3442 • Fax 301 695-7940 • www.dhgsonomed.com