

Hokanson MD35



- Bi-Directional Analog Doppler
 - Photo Plethysmography
 - Pneumo Plethysmography
- 2 speed Strip Chart Recorder
 - Arterial and Venous
 - Manometer
 - Standard Cuff Set
 - 5 Year Warranty
 - Economical

MD35 Tests & Applications	Doppler	Photo	Pneumo
Upper/Lower Extremity Spectral Doppler	X		
Upper/Lower Extremity Analog Doppler	X		
Ankle Brachial Index	X	X	X
Upper/Lower Extremity Segmental Blood Pressure	X	X	X
Blood Pressure in all Digits		X	
Digit and Limb Pulse Volume Recording			X
Penile Pressure and Plethysmography	X	X	X
Thoracic Outlet Syndrome		X	
Venous Capacitance			X
Maximum Venous Outflow			X
Venous Reflux		X	

Options Available

NIVP3 Software
 TD312 Cuff inflator
 MV10 Cuff Selector
 Additional Cuffs
 Carrying Case

Dimensions: Height 4.75" (12 cm) x Width 15.13" (38.4cm) x Depth 11.38" (28.9cm)



DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: DHGsonomed@aol.com
 800 669-3442 Fax 301 695-7940